

Arm and leg lift

HMI NO: 72792 | PRODUCT NO: BL-H



Our arm and leg lift, which we have great success with, is used in the handling and disinfection of arms and legs before operation, or other heavy lifting jobs on patients.

ARM AND LEG LIFT

Arm and leg lifter can be adjusted to various heights, so the user has the correct working position, which is important in relation to working environment. Very large lifting range of 75 cm. Easy to move around and with 4 brake wheels.

- Electrical height adjustment of Leg- and arm lift
- Powder-coated lift arm
- Powder coated frame
- Individual lockable wheels
- Battery solution with charger



TECHNICAL DATA

- Width 60 cm.
- Length 100 cm.
- Wheel size 75 mm.
- Elevation between 136-211 cm.
- Weight 31 kg.
- Max load 100 kg.
- Electrical voltage 230/24 V
- Density IP 50th
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk