



SPECIAL COUCHES

Designed to ensure optimum working conditions for the staff and maximum comfort for the user

We observe all legislation and government requirements, including the European Directive 93/42/EEC on medical devices. We are certified according to EN/ISO 13485, and all our products are CE-marked.

2-part ECCO/EKG examination couch

Couch with offset lifting system and opening in the lying surface to have good working condition for ECCO/ECG measurements. The opening on the lying surface is leaving space for electrodes. The lid for the opening can be parked, so that it doesn't interfere. The couch can be lowered for easy access from a wheelchair.



2-part ECCO/EKG examination couch

Product no.: **EKG-2000**

3-part endoscopy couch

Made with two back-rests, so the couch doesn't need to be unlocked and rotated when switching patients lying position.

Good hygienic conditions with lying surface cushions without seams and iron frames. Artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort.

The couch is with electric adjustment of height and both back rests. Trendelenburg and anti-Trendelenburg function.



3-part endoscopy couch

Product no.: **SØ40000**

3-part leg couch

With a fixed seat part, and the possibility to angle the back and leg part upward.

Lying surface with an artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort. There is chrome frame that protects cushions from collisions. The couch can be lowered for easy access from a wheelchair.



3-part endoscopy couch

Product no.: **U-2000**

HMI no.: **89897**

Shower couch, length 1900 or length 2400 with scissors or column lift

Care and washing in a shower couch with excellent working conditions. Bed guards can be positioned parallel to the bed surface and thereby facilitate the transfer.

Fixed or detachable headboards, mattress with good comfort and proper drainage with screwed on flex hose. The couch can be lowered for easy access from a wheelchair. Stable and with the best wheels that make it easy to move and with central brake system. Equipped with battery function so wiring can be avoided. Good hygiene mattress support formed in one piece light alloy.



Shower couch 1900, scissors lift / Column lift

Product no.: BRUSE-SAKS-1900 + BR-SØJLE-1900

HMI no.: 103220 + 108853

Shower couch 2140, scissors lift / Column lift

Product no.: BRUSE-SAKS-2140 + BR-SØJLE-2140

HMI no.: 103223 + 108920

2-part shower couch with lifted surface

Lifted lying surface with a tub under for collecting dirt and water. Bed guards. Good comfort and proper drainage with screwed on flex hose.

The couch can be lowered for easy access from a wheelchair. Stable and with the best wheels that make it easy to move around and with central brake system. Equipped with battery function so wiring can be avoided.



2-part shower couch with lifted surface

Product no.: BRUSE-SAKS-1000

HMI no.: 105502

Tilt couch with foot plate or divided adjustable footplates

Tilt function for rehabilitation, blood pressure tests and patients with syncope. Can be tilted 0-90 degrees and at 0 degrees function as a normal couch. Equipment rails on both sides. Degree indicator at 0 – 60 -90 degrees.

Emergency lowering function device can bring the couch surface to 0 degrees easily and quietly if everything fails. Can be equipped with armrests and straps. Central brake system.

Tilt couch with foot plate

Product no.: **VP1000-10**
HMI no.: **89899**

Tilt couch with adjustable footplates

Product no.: **VP1000**



Rehabilitation couch

Couch that can be tilted up to 45 degrees on the long side, so that the patient feels pressure under the feet. 2-part surface with back rest and mattress.

All four sides can be built up in two layers and are all build as cushions. The wheels are pulled up when the couch tilts so it stands stable and firm.

Rehabilitation couch

Product no.: **TR2-10000**
HMI no.: **79758**



rehab-care®

WE HELP YOU MAKE A DIFFERENCE



Rehab-Care develops, manufactures, services and repairs high-quality equipment and aids for the nursing and health sectors – e.g. hospitals, medical clinics and nursing homes.

At Rehab-Care the user's demands and needs are always in focus in the development and production of equipment and aids for hospitals, clinics and the nursing sector.

All Rehab-Care products are designed to ensure optimum working conditions for the staff and maximum comfort for the user.

We have detailed knowledge of all models, and we can customize our products to suit functional as well as aesthetic needs.

Our products need to be practical and functional, but we want them to look great too. If you have special wishes for e.g. upholstery colour, etc., we will customize it to suit your wishes and needs.

Your daily life should not be disrupted by long periods of interruption. We therefore service equipment on a day-to-day basis. We offer this service for both our own and third-party products.

All our prototypes are tested in co-operation with the Danish Technological Institute and other external consultants.

See all our products at rehabcare.dk