



ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	L-1000	Undervogn/saks	1
2	F7-VIP	Mellem ramme vip-leje	1
3	L8-A	Trd. travers spec	1
4	F1-VIP	Ben-del vip-leje	1
5	LA31 bund	Lift act	1
6	Act push rod	Dos not apply	2
7	F3-VIP	Ryg del	1
8	ISO 7380 - M10 x 16 --- 16N		4
9	F10	Udløser håndtag fys	1
10	F11	Forbindelsesled til Gas f.	1
11	F8	Holder for gasfjeder	1
12	GAS-02	For gasfjeder	2
13	GAS-17	Gasfjeder kort	1
14	F2	Hovedstykke	1
15	LA31-3 bund	LA31	1
16	CB12		1
17	F4-XX	Ryg plade	1
18	F5-XX	Nakkepude Fys-leje	1
19	U-51-polst-VIP		1
20	L38	Bremsebøjle	1
21	SHM5-VIP	Tilbehørs skinne, ryg	2
22	L28	Holder for gas fjeder	1
23	GAS-12-500N	Gasfjeder lang	1
24	L27	Udløserhåndtag	1
25	PAF 20165 P10	Permaglide	4
26	F-V-22	Led for leje	2
27	F-V-23		1
28	SHM3-VIP	Tilbehørsskinne	2
29	L53-H	Sliding block	2
30	L40	Plastbøsning	6
31	L31	Udløserhus trend	1
32	U20	Laske	1
33	L16-6	Montagebøsning	4
34	2474DYK125R26-28S45		4

Rev	Date	Note
		Scale: 1:20 Paper: A3 Sheet 2 of 2
Description: 3-delt FYS leje med vippe funktion		
Surface treatment:		Color:
Drawn By: MG	Drawn Date: 17-12-15	
Weight: 167438.70	State: - ok	Material:
 Tolerances ISO 2768-mK		Drawing No.: FYS-VIP Rev.