



ITEM NO.	PART NUMBER	Description	QTY.
1	SV-5000	Vogn for sengevipper	1
2	5945UAPT00P30-10 LH	Wheel low ø100 brake	4
3	L53-H	Glideklods	2
4	SV-1026	Trækstang	2
5	SV-6000	Løftearm	1
6	SV-1015	afstandsring	2
7	SV-1017	Krog for seng 5	2
8	Jumbo	CBJ styring	1
9	SV-7000	Skubbehåndtag	1
10	ISO 7380 - M6 x 10 - 10N		1
11	ZS-12x1-75	Pos bolt for adjust	2
12	DIN 912 M8 x 25 --- 25N		2
13	LA40-120	Linak LA40	1

Rev	Date	Note
A	04-03-17	

 Avnvej 10 - 7400-Herring Denmark WWW.rehab-care.dk	Scale:	1:5
	Paper:	A3
	Sheet 1 of 1	
Description: Senge vipper		
Surface treatment:	Males	Color: RAL 9006
Drawn By:	MG	Drawn Date: 03-09-18
Weight:	23602.57	State: - ok
		Material:
 Tolerances ISO 2768-mK	Drawing No.: SV-40000	Rev. A