



ITEM NO.	PART NUMBER	Description	QTY.
1	fro-001	Ryg PU pude	2
2	ts-6110	ryg beslag	4
3	ISO 7380 - M6 x 20 - 20N		8
4	ISO 7380 - M6 x 16 -- 16N		4
5	ts-7114	Holder til D-ryg	2
6	ISO 7380 - M6 x 16 - 16N		4
7	ISO 8678-M6x40-35-N		2
8	M6 greb	M6 greb	2
9	ts-7116	ryg plade	2
10	ts-7117	Samleplade	1
11	DIN 1587 - M6 --NUC		2

Rev	Date	Note

 rehab-care.dk[®] <small>Avnvej 10 -7400-Herring Denmark WWW.rehab-care.dk</small>	Scale: 1:2.5
	Paper: A3
	Sheet 1 of 1

Description: Ryg 2-delt PU komplet med beslag

Surface treatment:	Color:	
Drawn By: MG	Drawn Date: 07-03-19	
Weight: 5893.22	State: - ok	Material:

 Tolerances ISO 2768-cl	Drawing No.: TS-11300	Rev.
---	---------------------------------	------