




ITEM NO.	PART NUMBER	Description	QTY.
1	DRP-11000	Støttegreb venstre	1
2	DRP-12000	Støttegreb højre	1
3	Kugle ø50-M10 sort	Kugle ø50/M10	2



Rev	Date	Description
		Scale: 1:2.5
<small>Avnvej 10 -7400-Herring Denmark WWW.rehab-care.dk</small>		Paper: A3
Sheet 1 of 1		
Description: Oprejsning stag sæt		
Surface treatment:		Color:
Drawn By: MG	Drawn Date: 19-02-20	
Weight: 854.89	State: - ok	Material:
 Tolerances ISO 2768-cl	Drawing No.: DRP-22000	Rev. A