



ITEM NO.	PART NUMBER	Description	QTY.
1	ts-6050	Fodhviler bøjle	1
2	ts-7520	holder for knæpude	1
3	TST-2009	PU pude sidestøtte	1
4	ISO 7380 - M6 x 10 - 10N		6
5	ISO 7380 - M6 x 30 - 30N		1
6	DIN 1587 - M6 --NNU		1



Samles højre og venstre til et sæt

Rev	Date	Description	
 rehab-care.dk [®] Avnvej 10 -7400-Herning Denmark WWW.rehab-care.dk		Scale:	1:2.5
		Paper:	A3
		Sheet 1 of 1	
Description: knæpuder			
Surface treatment:		Color:	
Drawn By:	MG	Drawn Date:	25-03-2025
Weight:	1378.84	State: - ok	Material:
 Tolerances ISO 2768-cl		Drawing No.: ts-7530	Rev. A