

2-part column ECCO/EKG examination couch

PRODUCT NO: EKG-3000



2-PART COLUMN ECCO/EKG EXAMINATION COUCH

This examination couch on pillars is made with a customized lifting system and a smaller backrest in the lying surface to provide good working conditions when measuring EKG. The couch can be lowered significantly to make wheelchair access easier.

The backrest in the lying surface creates ample space for electrodes from below and can slide from side to side. The chrome frame protects cushions from impact. Skai withstands cleaning with alcohol and chlorine. The high-quality foam in the cushions ensures comfort.

- Electric height adjustment
- Very easy-to-clean wheel frame
- Powder-coated frame
- Lying surface made of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on a chrome tube frame to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels



TECHNICAL DATA

- Lying surface 200 x 70 cm
- Backrest height 46 cm
- Central brake wheels 125 mm
- Height of the couch between 65 and 105 cm
- Weight 80 kg
- Electric voltage 24 V
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk