

## 2-part examination couch, column lift

HMI NO: 86319 | PRODUCT NO: ULM1CBRA-SØ



Our good standard couch, is made to last. There are chrome frames that protects the cushions against collisions, artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort.

## 2-PART EXAMINATION COUCH, COLUMN LIFT

The couch can be lowered for easy access from a wheelchair. Large backrest angle adjustment for optimal seating position. Paper-holder and tear-off brace. Available with or without nose hole in cushion.

- Electric height setting
- Electric setting of Trendelenburg/antitrendelenburg
- Manual stepless position setting for back
- Very easy-to-clean wheel frame
- Comes with paper holder and tearing brace
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels

## TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Back section 80 cm
- Central brake wheels str. 125 mm
- Couch height between 50-90 cm (Can be changed to 65cm and 105cm)
- Weight 80 kg
- Maximum weight 250 kg
- Electrical voltage 24V
- CE marked

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

