

2-part examination couch, scissor lift

HMI NO: 72974 | PRODUCT NO: ULM1CBRA



Our good standard couch, is made to last. There are chrome frames that protects the cushions against collisions, artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort.

2-PART EXAMINATION COUCH, SCISSOR LIFT

The couch can be lowered for easy access from a wheelchair. Large backrest angle adjustment for optimal seating position. Paper-holder and tear-off brace. Available with or without nose hole in cushion.

- Electric height setting
- Manual stepless position setting for back
- Very easy-to-clean wheel frame
- Comes with paper holder and tearing brace
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Back section 80 cm
- Central brake wheels str. 125 mm
- couch height between 46-101 cm
- Weight 80 kg
- Maximum weight 250 kg
- Electrical voltage 24V
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

