

2-part gynecological couch

HMI NO: 89889 | PRODUCT NO: GYN-10000



Compact and stable 2-piece gynecological bed for leg holders. Does not require much space. Easy to access for patient and user. The wheel frame is made open at the front so that there is no frame in the way. The seat and back have stepless adjustment. Quality foam in the cushions ensures good comfort.

2-PART GYNECOLOGICAL COUCH

- Electric height adjustment
- Electric setting of the back
- Comes with paper holder
- Brackets for leg holders
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on a frame of chrome tubes to avoid damage to the cover
- 4 ø125 mm wheels with central brake

TECHNICAL DATA

- Cushion surface 120 x 70 cm
- Back section 90 cm
- Leg section 40 cm
- Central brake wheels ø125 mm
- Couch height between 48 and 90 cm
- Weight 50 kg
- Max. load 250 kg
- Electrical voltage 24 V
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

