

## 2-part rest- and transport bed, 2-layer

PRODUCT NO: AKUT-11000



2-piece rest and transport bed. A universal resting and transport bed with good comfort for patients and with focus on good ergonomics and stability for the staff. The bed can be driven far down, making access from a wheelchair easier.

Swinging bed rails can be laid down and swung under the lying surface. This allows the staff to get all the way to the edge of the mattress when treating the patient. Removable gables for transfer in all directions. Lying surface in compact laminate with mattress handlebars. The bed is with a double bottom, which is intended for X-ray.

## 2-PART REST- AND TRANSPORT BED, 2-LAYER

- Electrical height adjustment
- Electrical adjustment of the back position
- Manuel Trendelenburg function
- Easy to clean.
- Powder-coated
- Possible to purchase a loose mattress made of fire-retardant foam with a durable cover.
- Anti-slip on the back side.
- Equipped with fold-away bedrails.
- Central break on every wheel
- 2 directional wheels

## TECHNICAL DATA

- Mattress 200 x 70 cm
- Back support length 75 cm
- Back angle 85°
- Central break wheel size  $\varnothing$ 125 mm
- Bed height between 46 og 101 cm
- Weight 100 kg
- Max load 250 kg
- Electrical voltage 24 V
- CE-marked.

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

