

2-part rest- and transport bed on columns

PRODUCT NO: AKUT-12000 / 16000



2-piece rest and transport bed. A universal resting and transport bed with good comfort for patients and with a focus on good ergonomics and stability for the staff. The bed can be lowered far down, making access from a wheelchair easier.

Swinging bed rails, which can be laid down and swung under the lying surface. This allows the staff to get all the way to the edge of the mattress when treating the patient. Removable gables for transfer in all directions. Lying surface in compact laminate with mattress handlebars. This bed is also available with a width of 90cm. (akut-16000)

2-PART REST- AND TRANSPORT BED ON COLUMNS

- Electric height setting
- Electric back setting
- Electric setting of Trendelenburg/ antitrendelenburg
- Very easy-to-clean wheel frame
- Powder coated
- Possible to purchase a loose mattress made of fire-retardant foam with a durable cover.
- 2 directionally controlled wheels at the foot end
- Central brake on all wheels
- Fold-down bed rails with swing function
- Removable gables

TECHNICAL DATA

- Mattress 200 x 70 cm
- Back support length 75 cm
- Back angle 85°
- Central break wheel size \varnothing 125 mm
- Bed height between 50 og 90 cm (Can be changed to 65cm and 105cm)
- Weight 100 kg
- Max load 250 kg
- Electrical voltage 24 V
- CE-marked.

ADD-ONS

See a list of available spare parts at rehab-care.dk

