

2-part rest and transport bed with column lift, 2xbottom

PRODUCT NO: AKUT-13000



2-PART REST AND TRANSPORT BED WITH COLUMN LIFT, 2XBOTTOM

2-piece rest and transport bed. A universal rest bed and transport bed with good comfort for patients and with a focus on good ergonomics and stability for the staff. The bed can be driven far down, so access from a wheelchair will be easier.

Swing side rails, which can be laid down and swung under the lying surface. This allows the staff Unhindered access to the mattress edge when treating the patient. Removable gables for transfer in all directions. Lying surface in compact laminate with mattress handlebars. The bed is with a double bottom, which is intended for X-ray.

- Electric height adjustment
- Electric back adjustment
- Electric setting of Trendelenburg/antitrendelenburg
- Very easy to clean wheel frame
- Powder coated
- Possible to Purchase loose mattress of fire-retardant cold foam with durable cover. Anti-slip on the back.
- 2 guided directional wheels at the foot end
- Central brake on all wheels
- Fold-down bed rails with swing function
- Removable gables

TECHNICAL DATA

- Lying surface 200 x 70 cm
- Back part 75 cm
- Back angle 85°
- Central brake wheel ø125 mm
- Bed height between 50 and 90 cm (Can be changed to 65cm and 105cm)
- Weight 110 kg
- Max user weight 250 kg
- Electrical voltage 24 V
- CE-mark

ADD-ONS

See a list of available spare parts at rehab-care.dk



Avnvej 10, DK-7400 Herning | rehab-care.dk
Phone: +45 6915 3000 | E-mail: mg@rehab-care.dk