

## 2-part rest and transport bed with scissor lift

PRODUCT NO: AKUT-10000



2-part rest and transport bed. A universal rest bed and transport bed with good comfort for patients and with a focus on good ergonomics and stability for staff. The bed can be driven far down, so access from a wheelchair becomes easier. In addition, the bed is available with a width of 90 cm

## 2-PART REST AND TRANSPORT BED WITH SCISSOR LIFT

Swing-away rails can be laid down and swung under the sleeping surface. This allows the staff to get all the way to the edge of the mattress when treating the patient.

Removable gables for movement in all directions. Sleeping surface in compact laminate with mattress manage.

- Electric height adjustment
- Electric back adjustment
- Trendelenburg function
- Very easy to clean
- Powder coated
- Purchase a loose mattress made of fire-retardant cold foam with a durable cover. Anti-slip on the back.
- 2 pcs. directional wheels at the footing.
- Central brake on all wheels
- Fold-down bed rails with swing away function
- Removable gables



## TECHNICAL DATA

- Lying surface 200 x 70 cm
- Back part 75 cm
- Back angle 85°
- Central brake wheel ø125 mm
- Bed height between 46 and 101 cm
- Weight 100 kg
- Max user weight 250 kg
- Electrical voltage 24 V
- CE mark

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)



Avnvej 10, DK-7400 Herning | [rehab-care.dk](http://rehab-care.dk)  
Phone: +45 6915 3000 | E-mail: [mg@rehab-care.dk](mailto:mg@rehab-care.dk)