

2-part scan and x-ray couch on columns

PRODUCT NO: SØ-30000



Two-part scan bed / X-ray couch on column, which provides the best conditions for scanning tasks. The bed accommodates all conditions for staff and patients. Optimal space for scanning over the entire lying surface.

Lifting columns located at the outside. Longitudinal tube all the way to the floor. Very easy to access with bed rails folded under the lying surface. The

bed has central brakes, rounded edges and a flat surface in compact laminate. Removable gable hangers, back adjustment, trendelenburg function and fold-in bed rails.



2-PART SCAN AND X-RAY COUCH ON COLUMNS

- Electric height setting
- Electrical stepless adjustment of positions for back
- Electric setting of Trendelenburg/antiTrendelenburg
- Cobo control with removable jumbo battery
- Powder coated undercarriage
- Purchase option: mattress made of fire retardant cold foam with durable cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Madras surface 200 x 75 cm
- Back section 90 cm
- Back angle 60 °
- Central brake wheels size 125 mm
- Height between 50 and 90 cm
- Weight 110 kg
- Max Load 250 kg
- Electrical voltage 230/24v"
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

