

2-piece examination bed with scissor lift



The 2-piece examination bed with scissor lift from Rehab Care is made to withstand intensive daily use. Fitted with hydraulic pump that frees you from cables and wires. The chrome frame protects cushions from collision. The cover is made of artificial leather that can withstand cleaning with alcohol and chlorine. The good foam quality of the cushions ensures good comfort.

You can drive the bed so far down that access from a wheelchair becomes easier. There are great opportunities to adjust the back angle to give the patient optimal sitting position. Paper tear-off hanger helps hold the paper in place.

The couch is available with or without nose hole in the cushion and in different colors.

2-PIECE EXAMINATION BED WITH SCISSOR LIFT

- Hydraulic pump
- Manual stepless back position setting
- Highly easy to clean wheel frame
- Comes with paper holder and tear-off hanger
- Powder coated base
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on chrome tube frame to avoid damage to the cover
- 2 directional wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Back section 80 cm
- Central brake wheels str. 125 mm
- couch height between 46-101 cm
- Weight 80 kg
- Maximum weight 250 kg
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk