

3-part bariatric couch

PRODUCT NO: SØ80000



Couch with 2 back rests. Large flexible lying surfaces, as well as extra width. Four motors on the lifting function and a Trendelenburg function. Removable bed rails. Good hygiene conditions with lying surface cushions without stitching. Operated with hand control and/or foot control.

3-PART BARIATRIC COUCH

Good hygienic conditions with lying surface cushions without seams. Artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort. The couch is with electric adjustment of height and both back rests. Trendelenburg and anti Trendelenburg function.

- Electric height setting
- Electric adjustment of both backs
- Electric setting of Trendelenburg/antiTrendelenburg.
- Removable bed rails
- COBO system, battery function
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 210 x 90 cm
- Back portion 70 x 90 cm
- Back angle 85°
- Central brake wheels size 125 mm
- Height between 50 and 90 cm (Can be changed to 65cm and 105cm)
- Weight 150 kg
- Max Load 375 kg
- Electrical voltage 24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

