

## 3-part endoscopy couch

PRODUCT NO: SØ40000



Made with two back-rests, so the couch doesn't need to be unlocked and rotated when switching patients lying position.

## 3-PART ENDOSCOPY COUCH

Good hygienic conditions with lying surface cushions without seams and iron frames.

- Electric height setting
- Electric adjustment of both backs
- Electric setting of Trendelenburg/antiTrendelenburg
- Can be delivered with removable battery system
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- 2 directionally guided wheels
- Central brake on all wheels

## TECHNICAL DATA

- Cushion surface 210 x 70 cm
- Back portion 70 x 70 cm
- Back angle 85 °
- Central brake wheels size 125 mm
- Height between 50 and 90 cm (Can be changed to 65cm and 105cm)
- Weight 100 kg
- Max Load 250 kg
- Electrical voltage 230/24v
- CE – marked

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

