

3-part examination ml2



The 3-part examination bed has a lying surface with fixed gynecological seat and 2 identical adjustable backs.

The 3-part combi gynecological/ examination bed with scissor lift from Rehab Care is made to withstand intensive use every day. It is easy to clean the artificial leather cover, which can withstand cleaning with alcohol and chlorine. The good foam quality in the cushions ensures good comfort. The bed can be driven far down, so access from a wheelchair becomes easier. Large back angle adjustment for optimal sitting position.

Buy the bed with or without a nose hole in the back cushion. The seat is also available without a gyn cutout.

3-PART EXAMINATION ML2

- Electric height adjustment
- Manual stepless adjustment of positions for backs
- Very easy-to-clean wheel frame
- Powder-coated chassis
- Sleeping surface of fire-retardant cold foam with durable cover
- 4 pcs. wheels with brakes

TECHNICAL DATA

- Lying surface 200 x 70 cm
- Back part 80 cm
- Leg part 120 cm
- Brake wheels 125 mm
- Bed height between 46 and 101 cm
- Weight 80 kg
- Max load 250 kg
- Electrical voltage 24 V
- CE mark

ADD-ONS

See a list of available spare parts at rehab-care.dk

