

3-part multifunction couch on column

PRODUCT NO: ML1CBRA-SØ



3-PART MULTIFUNCTION COUCH ON COLUMN

The 3-part gynecological bed has a lying surface with gynecological seat and 2 identical backs. The couch is the most chosen bed among general practitioners. The 3-piece combi gynaecological/multifunction bed with column lift from Rehab Care is made to withstand intensive use every day. The chrome frame protects cushions from impact.

It is easy to clean the cover of artificial leather, which can withstand cleaning with alcohol and chlorine. The good foam quality of the cushions ensures good comfort.

The bed can be driven far down, so access from a wheelchair will be easier. Large back angle adjustment for optimal sitting position. Order the bed with or without nose hole in the back cushion.

- Electric height setting
- Electric setting of Trendelenburg/antitrendelenburg
- Manual stepless adjustment of positions for 2 backs and seat
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Lying surface 200 x 70 cm
- Back part 80 cm
- Leg part 120 cm
- Central brake wheels 125 mm
- Bed height between 50 and 90 cm (Can be changed to 65cm and 105cm)
- Weight 80 kg
- Max load 250 kg
- Electrical voltage 24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk



Avnvej 10, DK-7400 Herning | rehab-care.dk
Phone: +45 6915 3000 | E-mail: mg@rehab-care.dk