

3-part treatment couch, scissor lift

PRODUCT NO: FYS1CBRA



Our treatment couch is made to last. The chrome frames protects the cushions against impact, artificial leather in a quality that can withstand cleaning in alcohol and chlorine. The right foam quality cushions for best comfort.

The couch can be lowered for easy access from a wheelchair. Large

backrest angle adjustment for optimal seating position. High flexibility in the 3-part lying surface with 2-part back where conical head piece can be angled toward the floor. Available with or without nose hole.



3-PART TREATMENT COUCH, SCISSOR LIFT

- Electric height setting
- Manual stepless setting of positions for the head, back and legs
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Head part 35 cm
- Back section 50 cm
- leg part 120 cm
- Central brake wheels str. 125 mm
- Couch height between 46 and 101 cm
- weight 85 kg
- Maximum weight 250 kg
- Electrical voltage 24V
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

