

# 4-part multifunction couch on columns

PRODUCT NO: SOGYNLEJE



All function are possible with 4-part mattress support from gyn-seat, two-piece back, one with conical adjustable head piece .Our gynecological multi couch is made to last. There are chrome frame that protects the cushion against collisions, artificial leather in a quality that lasts for cleaning in alcohol and chlorine. The right foam quality cushions for comfort.

## 4-PART MULTIFUNCTION COUCH ON COLUMNS

The Couch can be lowered for easy access from a wheelchair. Large backrest angle adjustment for optimal seating position. Greater flexibility in the processing surface with 2-part back, where conical head piece can be angled toward the floor. Available with or without nose hole cushion.

- Electric height setting
- Electric setting of Trendelenburg and antitrendelenburg
- Manual stepless adjustment of positions for 2-piece main part, back and gyn seat
- Comes with paper holder
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels



## TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Main body 35 cm
- Upper back part 50 cm
- The middle back part 85 cm
- Gyn. The seat portion 30 cm
- Brake Wheel size 125 mm
- Couch Height between 50 and 90 cm (Can be changed to 65cm and 105cm)
- Weight 90 kg
- Max. Load 250 kg
- CE-mark
- Electrical voltage 24 V

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)