

4 part multifunction couch, scissors

PRODUCT NO: GYNLEJE



All functions are possible with 4-part cushions lying surface, gyn-seat prepared for leg support, two backs, one with conical adjustable head piece. Our gynecological multi couch is made to last. Chrome frames protects the cushions against collisions, artificial leather in a quality that can withstand cleaning in alcohol and chlorine. High quality foam in the cushions for best comfort

The Couch can be lowered for easy access from a wheelchair. Large backrest angle adjustment for optimal seating position. High flexibility in the 4-part lying surface with 2-part back, where conical head piece can be angled toward the floor. Available with or without nose hole in cushion.

Video 2 shows the bed with motor added in the back

4 PART MULTIFUNCTION COUCH, SCISSORS

- Electric height setting
- Manual stepless adjustment of positions for 2-piece main part, back and gyn seat
- Comes with paper holder
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on a frame of chrome pipe to avoid damage to the cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Main body 35 cm
- Upper back part 45 cm
- The middle back part 80 cm
- Gyn. The seat portion 30 cm
- Brake Wheel str. 125 mm
- Height between 46 and 101 cm
- Weight 90 kg
- Max. load 250 kg
- Electrical voltage 24 V
- CE-mark

ADD-ONS

See a list of available spare parts at rehab-care.dk

