

Amputations sling




The sling is used for lifting or moving a person with disability in hospitals, nursing homes, assisted living facilities, rehabilitation centers and in private.

AMPUTATIONS SLING

The sling is designed for both mobile hoists, ceiling hoist systems etc. and is ideally suited for lifting a person into and out of a seated position.

The sling is made for lifting an amputee person with reduced function in the upper body, and femoral hip joint, but with control over the head. The sling is floor reaching

TECHNICAL DATA

- Material: Polyester
- Max load: 275 kg.
- Washing: 

ADD-ONS

See a list of available spare parts at rehab-care.dk

