

# Bed carrier



The Bed Carrier facilitates the transfer of beds and to promote good working environment. Used in nursing homes and in institutions.

## BED CARRIER

The Bed Carrier is motor-driven, has a removable battery and is provided with four locking wheels.

- Electric tilt
- Lockable wheels
- Battery powered

## TECHNICAL DATA

- Length 120 cm
- Width 82 x 132 cm
- Angle 60°
- Min free space under bed 155 mm
- Max load 150 kg
- Wheel diameter 100 mm
- waterproof class IPX4
- CE-marked

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

