


Bed lift (Wide)

PRODUCT NO: SE-140-200



BED LIFT (WIDE)

The bed lift is to be used for the consumers own private bed. Therefore, it can be adjusted for all types of elevation beds. The bed lift can be mounted under the bed, and the lifting system is known from the hospitality beds.

High lifting height for correct working position. Stable lift from the bottom and up. Removes the need of a nursing bed in the home, and it can be mounted on one bed so that the double bed can be preserved.

- The bed lift can be mounted on the bed or it can be delivered as a complete bed with a frame.
- Adjustable for widths between 120-140cm
- Stable scissor lift
- Electrical height adjustment
- Lockable wheels

TECHNICAL DATA

- Length 120cm
- Minimum bed width 120cm
- Maximum bed width 140cm
- Height min 17cm
- Height max 75cm
- Lifting interval 58cm
- Maximum load 200kg
- User weight 170kg
- Weight 40kg
- Wheel diameter 100 mm
- Density class for water 66 IP

ADD-ONS

See a list of available spare parts at rehab-care.dk

