

# Bed turner

PRODUCT NO: SV400



Hospital bed service trolley for servicing, maintenance and cleaning. Lifts all bed types from the floor and can further tilt the bed so, that all work on the bed may be done in correct working position.

## BED TURNER

Easy to operate and move round. The levers can be adjusted in all directions and turned, so lift underneath the lying surface also is possible.

- Electrical height adjustment
- Powder-coated undercarriage
- Brake function on wheels
- Adjustable to bed format width and depth



## TECHNICAL DATA

- Width 69-100 cm.
- Height 44-60 cm.
- Weight 80 kg.
- Max. Load 220 kg.
- height between 63 – 212 cm.
- IP class IP X6
- CE-Marked

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

