

Bobath couch

HMI NO: 79750 | PRODUCT NO: B100-70



Rehabilitation couch that has the best possible properties and durability in daily use. The couch can be lowered for easy access from a wheelchair. Artificial leather in a quality that can withstand cleaning in alcohol and chlorine. Stable and with the best wheels that make it easy to move and with central brake. The lying surface can be made in the desired dimensions of the customer.

BOBATH COUCH

- Electric height setting
- Very easy-to-clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- 2 directionally guided wheels
- Central brake on all wheels

TECHNICAL DATA

- Cushion surface 200 x 70 cm
- Height between 38 og 92 cm
- Lifting interval 54 cm
- Brake wheel 125 mm
- Weight 75 kg
- Max. load 250 kg
- Electrical voltage 24 V
- IP klasse IPX4
- CE-Marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

