

Chair for blood sampling

PRODUCT NO: DRP-3000



A functional chair for blood sampling, developed together with bioanalysts. The column lift allows staff to perform blood collection sitting or standing. Ergonomically correct solution.

The structure of the chair is designed for many daily blood collections without wear and tear. Equipped with removable battery. The seat part can be rotated 120 degrees relative to the wheel frame. The chair is available with and without central brake

CHAIR FOR BLOOD SAMPLING

- Electric height adjustment
- Battery function
- Possibility of central braking system

TECHNICAL DATA

- Length 56 cm
- width 48 cm
- Brake wheel 100 mm
- Height between 50 – 90 cm
- Weight 35 kg
- Max. load 200 kg
- Electrical voltage 24V
- CE-mark

ADD-ONS

See a list of available spare parts at rehab-care.dk

