

Comfort sling with head rest, mesh

HMI NO: 115415 | PRODUCT NO: 40049



The sling is used for lifting or moving a person with disability in hospitals, nursing homes, assisted living facilities, rehabilitation centers and in private.


COMFORT SLING WITH HEAD REST, MESH

With the open form of the sling it is possible to place the sling under the patient's buttocks, making the installation and removal comfortable for the patient and easy for the caregiver.

The sling is designed for both mobile hoists, ceiling hoist systems etc. and is ideally suited for lifting a person into and out of a seated position.

The sling is made for lifting a person with reduced function in the upper body, and femoral hip joint, but with control over the head. The sling is floor reaching.

TECHNICAL DATA

- Material: Polyester
- Max load: 275 kg.
- Washing: 

ADD-ONS

See a list of available spare parts at rehab-care.dk

