

DayCare chair

HMI NO: 72976 | PRODUCT NO: DSL-20000 + DSLHYNDEGRAA



DayCare chair designed for day surgery departments and dialysis. Used as therapy chair, examination and treatment chair.

Best seating comfort and many adjustment possibilities. Good handling, with fold-away and removable armrests. Everything is made resilient for daily use and can easily be kept clean. Developed in close collaboration with staff in hospitals.

DAYCARE CHAIR

We are very responsive, if there are good ideas for improvements. The DayCare chair is used by hospitals and the healthcare sector. Used as a transport chair between surgery and recovery or for prior study of patients. The couch can be equipped with many accessories as it is widely used in many departments.

- Electrical adjustment of all functions of the chair, such as lift, Trendelenburg, back, leg, footrest part
- Powder-coated undercarriage
- Removable armrest, controlled parallel together with seat part
- Surface cushion of anti-inflammable foam covered with hard-wearing material

TECHNICAL DATA

- External dimensions with armrests 235 x 89 cm
- Cushion surface 210 x 60 cm
- Central brake wheels str. 125 mm
- Seat height from 57 to 93 cm
- Weight 95 kg
- Max Load 200 kg
- Electrical voltage 24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

