

Emmy patient lift

HMI NO: 115025 | PRODUCT NO: PS-46000



Emmy patient lift, vertical lift with 3-part column, is for the stable and comfortable transfers. Equipped with the latest technology and the best components, to ensure the best durability and function.

EMMY PATIENT LIFT

A hoist which is easy to work with, since it has the best quality wheels and parts made for good handling.

- Electrical adjustment by remote control.
- Electrical angling of legs.
- Emergency lowering.
- Emergency stop.
- Removable battery.

TECHNICAL DATA

- Width 70 cm.
- Length 118 cm.
- Height 150 cm.
- Weight 45 kg.
- Max. Load 150 kg.
- Height between 26 and 176 cm.
- IP code IPX4
- CE marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

