

EPJ table

HMI NO: 100712 | PRODUCT NO: DRP-650



The EPJ tables are simple, but also very stable as they are made 10 mm, stainless steel plate. This makes the tables heavier so that they do not inadvertently fall over.

Very easy operation with up-down control panel, and a removable battery. There are also cutouts in the table top for any power cable. Electric height adjustment

EPJ TABLE

- Brake function on all wheels
- Handle in the table for easy placement
- Removable battery

TECHNICAL DATA

- Table surface 55 x 45 cm
- Height between 65 – 132 cm
- Brake wheels 100 mm
- Weight 25 kg
- Max. Load 50 kg
- Power 230/24 V
- CE – mærket

ADD-ONS

See a list of available spare parts at rehab-care.dk

