

ErGo-Basi 1000

HMI NO: 108699 | PRODUCT NO: TS-40000



Toilet- and showerchair with manual height adjustment.

Simple and functional design. The toilet- and shower chairs from Rehab-Care are hygiene optimized, with a minimum of horizontal surfaces and with a limited time-consumption for cleaning. The chair is easily run and is occupying only a limited space. Made from good materials for long life span. Customizable. Large accessory catalog

ERGO-BASI 1000

- Soft PU seat with integrated bedpan rails.
- Seat height manual adjustable, in 5 steps.
- Armrests height adjustable and fold-up.
- Back in Skai with Velcro, depth adjustable.
- Footrests strong and stable, height adjustable with fold-up footplates.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.
- Development and production in Denmark, subject to certified ISO13485.

TECHNICAL DATA

- Length 950 mm
- Width 580 mm
- Seat height 475 – 600 mm
- Distance between armrests 473 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle 0 deg.
- Castors with brakes, size Ø125 mm
- Weight of the chair 15 kg.
- Max load 200 kg.
- Electrical voltage 0V
- ISO17966:2016 certified

*Complete technical data sheet on request.

ADD-ONS

See a list of available spare parts at rehab-care.dk

