

ErGo OpTi

HMI NO: 134843 | PRODUCT NO: TS-37500



For dressing and undressing of the user. Hygiene optimized, with a minimum of horizontal surfaces, easy to keep clean. The toilet-bath chair is easy to use and takes up very little space. Optimal working conditions for helpers. Ergonomic comfort for the user. Made of good materials for the best service life. Custom options for specific users. Large catalog of accessories.

ERGO OPTI

- Flexible settings ensure optimal working position.
- Dressing and undressing the user
- Seat height electrically adjustable.
- Seat tilt electric.
- Armrests foldable.
- Head support height and depth adjustable.
- Foot rests strong and stable, height adjustable with foldable footplates.
- Battery indicator in hand control.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.
- Development and production in Denmark, subject to certified ISO13485.

TECHNICAL DATA

- Seat height 51-91cm
- Seat angle 0-90 degrees
- Seat width 43cm
- Seat depth 42cm
- Toilet seat opening 18cm
- Hygiene opening 12cm
- Total width 58cm
- Total depth 100cm
- Height in lowest position 10cm
- Weight og the chair 35 kg
- Castors with brakes, size Ø125 mm
- Max load 200 kg
- Electrical voltage 24v
- CE marked
- IPx6W
- ISO17966:2016 certified

ADD-ONS

See a list of available spare parts at rehab-care.dk

