

# ErGo-SelPro 1500

HMI NO: 111000 | PRODUCT NO: TS-39700



**Toilet- and shower chair with manual height adjustment. 2 X 24" drive wheels and 2 X 125 mm frontcastors.**

The toilet- and shower chairs from Rehab-Care are hygiene optimized, with a minimum of horizontal surfaces and with a limited time-consumption for cleaning. The chair is easily run and is occupying only a limited space.

## ERGO-SELPRO 1500

- Soft PU-seat with integrated bedpan rails.
- Armrests adjustable in height and foldable.
- Backrest in Skai with Velcro, adjustable in depth.
- Footrests strong and stable, height adjustable with foldable footplates.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.

## TECHNICAL DATA

- Length 1014 mm
- Width 576 mm
- Seat height mellem 460 – 520 mm
- Distance between armrests 473 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle 0 gr./0 gr.
- Castors with brakes, size Ø125 mm or 2 x 24" wheels
- Weight of the chair 20 kg
- Max load 200 kg
- Electrical voltage 0V
- ISO17966:2016 certified

\*Complete technical data sheet on request.

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)



®