

ErGo-Tilt 3000

HMI NO: 107344 | PRODUCT NO: TS-38000



Toilet- and shower chair, manual height adjustment, electrical tilt-function.

The toilet- and shower chairs from RehabCare are hygiene optimized, with a minimum of horizontal surfaces and with a limited time-consumption for cleaning. The chair is easily run and is occupying only a limited space.

ERGO-TILT 3000

- Soft PU-seat with integrated bedpan rails.
- Armrests adjustable in height and foldable.
- Backrest in Skai with Velcro, adjustable in depth.
- Footrests strong and stable, height adjustable with foldable footplates.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.

TECHNICAL DATA

- Length 950 mm
- Width 580 mm
- Seat height 475 – 600 mm
- Distance between armrests 473 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle -5 gr./25 gr.
- Castors with brakes, size Ø125 mm
- Weight of the chair 17 kg
- Max load 200 kg
- Electrical voltage 24 V
- IP x6w
- CE marked
- ISO17966:2016 certified

*Complete technical data sheet on request.

ADD-ONS

See a list of available spare parts at rehab-care.dk

