

# ErGo-Tilt 4000

HMI NO: 108696 | PRODUCT NO: TS-25000



**Toilet- and showerchair, electrical height adjustment, electrical tilt-function.**

The toilet- and shower chairs from RehabCare are hygiene optimized, with a minimum of horizontal surfaces and with a limited time-consumption for cleaning. The chair is easily run and is occupying only a limited space.

## ERGO-TILT 4000

- Soft PU-seat with integrated bedpan rails.
- Armrests adjustable in height and foldable.
- Backrest in Skai with Velcro, adjustable in depth.
- Footrests strong and stable, height adjustable with foldable footplates.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.

## TECHNICAL DATA

- Length 976 mm
- Width 589 mm
- Seat height 520 – 900 mm
- Distance between armrests 473 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle -5° / +30°
- Castors with brakes, size Ø125 mm
- Weight of the chair 30 kg
- Max load 200 kg
- Electrical voltage 24 V
- IPx6w
- CE marked
- ISO17966:2016 certified

\*Complete technical data sheet on request.

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

