

Ergo tilt 4000 hydraulisk

HMI NO: 138661 | PRODUCT NO: TS-26500



Toilet-shower chair with hydraulic height adjustment and manual tilt via double-acting gas spring.

Simple and functional design. Hygiene optimized, with a minimum of horizontal surfaces, easy to keep



Avnvej 10, DK-7400 Herning | rehab-care.dk
Phone: +45 6915 3000 | E-mail: mg@rehab-care.dk

clean. The chair is equipped with hydraulic pump, is easily driven and takes up very little space. Optimal working conditions for the caregiver. Ergonomic comfort for the user. Made of good materials for best longevity. Customizable. Large accessories catalog.



ERGO TILT 4000 HYDRAULISK

- Seat in soft PU with integrated pelvic sliders.
- Seat height hydraulically adjustable.
- Seat tilt via gas spring, operated by grip at the top of the back.
- Armrests height adjustable and foldable.
- Back in Skai with Velcro, depth adjustable.
- Headrest height and depth adjustable.
- Footrests strong and stable, height adjustable with folding footplates.
- Withstands cleaning with rubbing alcohol and chlorine.
- Machine washable.
- Development and production in Denmark, subject to certified ISO13485

TECHNICAL DATA

- Length 976 mm
- Width 589 mm
- Seat height 520 – 900 mm
- Distance between armrests 473 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle -5° / +30°
- Castors with brakes, size Ø125 mm
- Weight of the chair 30 kg
- Max load 200 kg
- CE marked
- ISO17966:2016 certified

*Complete technical data sheet on request.

ADD-ONS

See a list of available spare parts at rehab-care.dk





Avnvej 10, DK-7400 Herning | rehab-care.dk
Phone: +45 6915 3000 | E-mail: mg@rehab-care.dk