

ErGo-Tilt 6000

HMI NO: 115024 | PRODUCT NO: TS-25300



Toilet- and showerchair, electrical height adjustment, electrical tilt-function.

The toilet- and shower chairs from RehabCare are hygiene optimized, with a minimum of horizontal surfaces and with a limited time-consumption for cleaning. The chair is easily run and is occupying only a limited space.

ERGO-TILT 6000

- Soft PU-seat with integrated bedpan rails.
- Armrests adjustable in height and foldable.
- Backrest in Skai with Velcro, adjustable in depth.
- Footrests strong and stable, height adjustable with foldable footplates.
- Withstands cleaning with alcohol and chlorine.
- Machine washable.

TECHNICAL DATA

- Length 843 mm
- Width 589 mm
- Seat height 520 – 900 mm
- Distance between armrests 353 mm
- Seat width 430 mm
- Seat depth 430 mm
- Seat angle -5° / +30°
- Castors with brakes, size Ø125 mm
- Weight of the chair 30 kg
- Max load 200 kg
- Electrical voltage 24 V
- IPx9
- CE marked
- ISO17966:2016 certified

*Complete technical data sheet on request.

ADD-ONS

See a list of available spare parts at rehab-care.dk

