

# Flexi Mover Base



## FLEXI MOVER BASE

The stand-up and transfer aid Mover is a great alternative to the wheelchair, as it offers comfortable sitting due to its padded swivel seat plates and knee supports. The Mover allows quick, stable and safe transport of the patient from the bed to a chair or toilet or vice versa. The rear castors can be fixed by means of the brakes. The Mover is useful for the transportation of patients, who have difficulties walking or standing.

## TECHNICAL DATA

### ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)