

Full sling







The sling is used for lifting or moving a person with functional reduction in hospitals, nursing homes, assisted living facilities, rehabilitation centers, riding schools and in private.

FULL SLING

RehabCare full sling is particularly suitable for people who are unstable in the neck area. The full sling can be placed under the armpits or support the body completely. A safe lifting sling also for leg amputated because the straps can be crossed. The sling is floor reaching.

TECHNICAL DATA

- Max. load: 275 kg.
- Material: polyester
- Washing:     

ADD-ONS

See a list of available spare parts at rehab-care.dk

