

FYS-hammel

HMI NO: 134845 | PRODUCT NO: FYS-HAMMEL



FYS-HAMMEL

Rehabilitation bed 70cm wide

The rehabilitation bed with scissor lift from Rehab Care is made to withstand intensive daily use. The chrome frame protects cushions from impact. The cover is made of artificial leather, which can withstand cleaning with alcohol and chlorine. The good foam quality of the cushions ensures good comfort.

Access from a wheelchair becomes easier when the bed is at its lowest point. There are great adjustment opportunities to the angle of both the main part and the leg part to give the patient optimal sitting position for rehabilitation.

- Electric height adjustment
- Electrical setting of back/head position
- 2 Actuators with mechanical emergency pull
- Manual stepless setting of position for legs
- Highly easy to clean wheel frame
- Powder coated undercarriage
- Lying surface of fire-retardant cold foam with durable cover
- Lying surfaces are mounted on chrome tube frames to avoid damage to the cover
- 2 directional wheels
- Central brake on all wheels

TECHNICAL DATA

- Lying surface 200 x 70 cm
- Back part 100 x 70
- Leg 45 x 70
- Seat part 45 x 70
- Bed height between 47cm – 101cm
- Max load 250kg
- CE marked
- Central brake wheel 125mm
- Weight 80kg
- Electrical voltage 24v

ADD-ONS

See a list of available spare parts at rehab-care.dk

