

# Gaja chair

PRODUCT NO: DRP-29000



Gaja Treatment chair for hospitals and healthcare sector. The chair is equipped with central lifting that does not take up unnecessary space around the user. Adjusting the height, back, and legs is with wireless bluetooth hand control or foot control, so wires do not get in the way. Different armrests can be selected to perfectly suit the user's needs in everyday life.

## GAJA CHAIR

- Flexible settings ensure optimal working position.
- Centrally located column motor
- BlueTooth controls
- Electric setting, up/down and tilt
- Motor functions for all lying surfaces
- Various armrest options. Configurable in many ways

## TECHNICAL DATA

- Seat height 55-88cm
- Seat width 60cm
- Total width 81cm
- Max load 200kg
- Trendelenburg 0° / 12°
- Seat depth 52cm
- Rotatable 180 degrees

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

