


Hoist unit AR-200

HMI NO: 100198



AR-200 hoist unit is for the smart, easy and silent transfer of patients. also presented in detachable version.

HOIST UNIT AR-200

TECHNICAL DATA

- Function: Patient moving and handling
- Capacity: (swl) 200 kg
- Dimensions: 306x270x148 mm
- Weight: 10 kg
- Power: Battery charged, 24 V 2,9 Ah
- Lifting Tape: 35x2400 mm
- Lifting velocity: 35-70 mm/sec
- Lifts per charge: >40 @120 kg
- Charging: By handset on magnetic wall-charger
- Charging time: 8 Hours
- Unique features:
 - Advanced battery protection system
 - Auto stand-by
 - No power consumption in stand-by

ADD-ONS

See a list of available spare parts at rehab-care.dk

