


Hoist unit AR-260

HMI NO: 100199



AR-260 hoist unit is for the smart, easy and silent transfer of patients. a light weight unit with a capacity for heavier loads.

HOIST UNIT AR-260

TECHNICAL DATA

- Function: Patient moving and handling, heavy duty
- Capacity (swl): 260 kg
- Dimensions: 377x298x135 mm
- Weight: 14 kg
- Power: Battery charged, 24 V 5 Ah
- Lifting Tape: 50x2400 mm
- Lifting velocity: 35-70 mm/sec
- Lifts per charge: >50 @ 150 kg
- Charging: By handset on magnetic wall-charger
- Charging time: 8 Hours
- Unique features:
 - Advanced battery protection system
 - Auto stand-by
 - No power consumption in stand-by

ADD-ONS

See a list of available spare parts at rehab-care.dk

