

# Hoist unit AR-500

HMI NO: 100200



AR-500 hoist unit is for the smart, easy and silent transfer of patients. The AR-500 unit is build for Bariatric patients and can take loads up to 500 kg.

## HOIST UNIT AR-500

## TECHNICAL DATA

- Function: Bariatric Patient moving and handling
- Capacity: (swl) 500 kg
- Dimensions: 495x337x188 mm
- Weight: 22 kg
- Power: Battery charged, 24 V 10 Ah
- Lifting Tape: 50x2400 mm
- Lifting velocity: 50-70 mm/sec
- Lifts per charge: >60 @500 kg
- Charging: By handset on magnetic wall-charger
- Charging time: 8-14 Hours
- Unique features:
  - Advanced battery protection system
  - Auto stand-by
  - No power consumption in stand-by

## ADD-ONS

See a list of available spare parts at [rehab-care.dk](http://rehab-care.dk)

