

Hospital table / Back table

PRODUCT NO: DRP-21000



The hospital table from Rehab-Care is available in different sizes. It is very easy to clean. The table is equipped with good wheels with brakes. Very easy to control with an up-down button.

- Can be used as a back table during surgery
- Electrical height adjustment
- Removeable battery
- Break function on all wheels

Hospitalsbord 1200 x 750. Varenr. **DRP-21000**

Bagbord 1000 x 600. Varenr. **DRP-20000**

HOSPITAL TABLE / BACK TABLE

TECHNICAL DATA

- Table surface 120 x 75/100x60cm or special dimensions
- Height between 71 – 112 cm
- Brake wheels 125 mm
- Weight 37 kg
- Max. Load 100 kg
- Electrical voltage 24V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

