

Hospital table

HMI NO: 89926



The hospital table from Rehab-Care is available in different sizes. It is very easy to clean. The table is equipped with good wheels with brakes. Very easy to control with an up-down button.

HOSPITAL TABLE

- Can be used as a back table during surgery
- Electrical height adjustment
- Removeable battery
- Break function on all wheels

Models

Hospital table 550 x 750. Product no.
DRP-710

Hospital table 1200 x 750. Product no.
DRP-21000

Hospital table 1000 x 600. Product no.
DRP-20000

TECHNICAL DATA

- Table surface 55 x 75 cm
- Height between 71 – 112 cm
- Brake wheels 100 mm
- Weight 37 kg
- Max. Load 100 kg
- Power 24 V
- CE – marked

ADD-ONS

See a list of available spare parts at rehab-care.dk

